



**PRODUCT
CATALOGUE**



ABOUT US



Chiron Surgimed LLP is a pioneering company at the forefront of innovative medical technology. Established in 2013, we have been dedicated to revolutionizing healthcare by developing cutting-edge medical devices that enhance patient care and improve clinical outcomes.

Our company specializes in the design, manufacturing, and distribution of a wide range of medical devices catering to various medical specialties such as Infusion Therapy, Respiratory Care, Urology and more.

In essence, **Chiron Surgimed LLP** is not just a company; it's a commitment to improving lives through innovative medical technology.

MEET THE TOP BRASS



Mr. Naresh Patel

Founder & Chairman

Founder to several companies with background in speciality polymer and in-depth experience of over 3 decades. His core competency is plastic polymers in the field of medical devices. Extremely diligent and hardworking, he is highly devoted spiritually which reflects in his day-to-day activities. Believes in working hard and simple living.



Mr. Abhay Kumar

CEO & Managing Director

He possesses strong leadership qualities with over 35 years of professional experience in leading global medical device companies. He is a post graduate in Microbiology and done business management from premiere institute. He has strong knowledge of global business development, product development, business management and quality management skills.



Mr. Harsh Patel

Director

A young entrepreneur with zeal to excel in Business Development and Marketing. Graduate in Commerce and keenly involved in Business Development activities and Production. Extremely hard working and intelligent person who wants to take Chiron to the next level. Talented and passionate about learning all the tricks of the trade.





INFUSION THERAPY

Infusion therapy is when you receive medication through a needle or catheter, usually intravenously (IV).

CHIRO-DRIP

IV Set Vented



Air Vented Spike



Roller Clamp



Male Luer
with Hypodermic Needle



chironsurgimed.com




- › Sharp air vented spike.
- › Suitable for pressure & gravity infusions.
- › Latex-free.
- › DEHP-free.

CHIRO-DRIP

IV Set Non-Vented



Non-Winged




Winged



CHIRO-DRIP Y

IV Set Vented with Y Injection Site



Y Injection Site



Air Vented Spike




- › Sharp air vented spike.
- › Suitable for pressure & gravity infusions.
- › Latex-free.
- › DEHP-free.

CHIRO-DRIP O

Oncology Infusion Set



Y Injection Site



Air Vented Spike



› Oncology Infusion Set are designed for the administration of chemotherapy.



CHIRO-DRIP MICRO

Micro IV Set Non- Vented



Micro Drip



Roller Clamp



Male Luer
with Hypodermic Needle



chironsurgimed.com



- › Micro I.V. Infusion set for paediatric use
- › Approximately 60 drops/ml
- › Sharp Spike for easy insertion

CHIRO-DRIP MICRO Y

Micro IV Set Vented with Y Injection Site



Y Injection Site



Micro Drip



chironsurgimed.com

- › Micro I.V. Infusion set for paediatric use
- › Approximately 60 drops/ml
- › Y-injection port designed for multiple injections
- › Sharp and Winged Spike for easy insertion

CHIRO-GUARD

IV Set Vented with 0.2 Micron Filter with Y Injection Site




0.2 Micron Filter



Roller Clamp



chironsurgimed.com

- 
- › Sharp air vented spike
 - › Suitable for pressure & gravity infusions
 - › Latex free
 - › DEHP free
 - › Y-injection port designed for multiple injections
 - › 0.2 Micron Inline Filter



CHIRO-FLOW

Flow Regulator



Flow Regulator intergrated with IV Set



Flow Regulator

- › Ergonomic dial design.
- › Compatible with standard sets.
- › Y-injection port designed for multiple injections.
- › DEHP-Free soft and kink resistant PVC tubing.

CHIRO-BT

Blood Transfusion Set



Non-Vented



Vented

CHIRO-VOL

Measured Volume Burette Set



Air Vented Spike



Silicone Floater



- › Burette chamber of 110ml/150ml.
- › Micro drip with drop size of 60 drops per ml.
- › DEHP-Free soft and kink resistant PVC tubing.



CHIRO-VOL Y

Measured Volume Burette Set with Y Injection Site



Y Injection Site



Silicone Floater




- › Burette chamber of 110ml/150ml.
- › Micro drip with drop size of 60 drops per ml.
- › Y-injection port designed for multiple injections.
- › DEHP-Free soft and kink resistant PVC tubing.



CHIRO-VEIN AND NEO

IV Cannula with Wings and Injection Port




TECHNICAL SPECIFICATION

Size	Catheter OD(mm)	Catheter length(mm)	Flow Rate (Water) (ml/min)	Box Qty. (Inner/Outer)
14G	2.1	45	270	100/1000
16G	1.7	45	180	100/1000
17G	1.5	45	125	100/1000
18G	1.3	45	80	100/1000
20G	1.1	32	54	100/1000
22G	0.9	25	33	100/1000
24G	0.7	19	20	100/1000
26G	0.6	19	19	100/1000

CHIRO-VEIN SAFETY

Safety IV Cannula with Wings and Injection Port




TECHNICAL SPECIFICATION

Size	Catheter OD(mm)	Catheter length(mm)	Flow Rate (Water) (m/min)	Box Qty. (Inner/Outer)
18G	1.3	45	80	50/500
20G	1.1	32	54	50/500
22G	0.9	25	33	50/500
24G	0.7	19	20	50/500

CHIRO-VEIN NEO SMALL WINGS


IV Cannula with Small Wings and without Injection Port



TECHNICAL SPECIFICATION

Size	Catheter OD(mm)	Catheter length(mm)	Flow Rate (Water) (ml/min)	Box Qty. (Inner/Outer)
24G	0.7	19	20	100/1000
26G	0.6	19	19	100/1000

3-WAY STOPCOCK



TECHNICAL SPECIFICATION


Options			Packaging Specification
Tap Colour	Material	Tap Rotation	inner
Red	Lipid Lipid Resistant	Smooth 45° Click	50
White	Lipid Lipid Resistant	Smooth 45° Click	outer
Blue	Lipid Lipid Resistant	Smooth 45° Click	1000

EXTENSION TUBES



**Extension Tube
with Filter**


**3 Way Stopcock
with Extension Tube**



Spiral Extension Tube



EXTENSION TUBES




Standard Extension Tube



Available in single, double and triple configuration.

Multi-Infusion Extension Tube with Needle Free





UROLOGY & NEPHROLOGY


A nephrologist deals with problems of the kidney with medicines and a urologist deals with kidney ailments that require surgery. Any ailments in the urinary system, prostate, bladder and penis which require surgical intervention are tended to by a urologist.

CHIROCATH

Foley Catheter



Latex Foley Catheter



- Sizes Available :**
- 8FR
 - 10FR
 - 12FR
 - 14FR
 - 16FR
 - 18FR
 - 20FR
 - 22FR



Silicon Foley Catheter



- Sizes Available :**
- 6FR
 - 8FR
 - 10FR
 - 12FR
 - 14FR
 - 16FR
 - 18FR

CHIRO-BAG

Urine Bag



Premium version features a high-quality plastic handle



- › Urine collecting bag used for urine drainage.
- › Two litre capacity bag provided with super smooth kink resistant tubing.
- › Efficient non-return valve with top outlet.

chironsurgimed.com



CHIRO-METER

Urometer




chironsurgimed.com



- › Graduated measured volume meter for accurate measurement of the urine output.
- › Urometer is directly attached to a 2 litre capacity urine bag having overflow facility.
- › Soft and kink resistant PVC tubing.



URETHRAL TUBE




Sizes Available :
from 8FR to 16FR



- > Designed for short term bladder catheterization through urethra.
- > Manufactured from medical grade PVC compatible with catheter lubricants.



DJ STENT



- > Short and Long Term connection of Urethra and Bladder.
- > Closed Tip and Open Tip

HAEMODIALYSIS KIT

Haemodialysis Double/Triple Lumen Catheter



chironsurgimed.com

› Haemodialysis catheter are double/ triple lumen catheter that provide temporary vascular access for haemodialysis until a permanent access is available or until another type of dialysis therapy is substituted.

› **Kit Contents :**

- Indwelling catheter
- Y-Introducer needle with valve
- J-tip guide wire
- Dilators- 2 Nos
- Luer lock syringe
- Scalpel



RESPIRATORY CARE

Respiratory therapy is a field of treatment that helps individuals optimize their respiratory function, breathe more easily, and live more comfortably.

AMBU BAGS

Silicone Resuscitator (Adult/Pediatric/Neonatal)



- > Silicone resuscitator (except mask, oxygen tubing and reservoir bag) can be autoclaved repeatedly.
- > It is semi-transparent and come with a pressure limitation valve for patient safety.

GUEDEL AIRWAY




chironsurgimed.com



- > The center passage permits the passage of suction catheters to clear secretions from the back of the throat.




Sizes Available :



CHIRO-OXY


Oxygen Mask (Adult/Pediatric/Neonatal)



- › The clinically clean oxygen masks are used to administer oxygen.
- › The adjustable nose clip and the under-chin design increase the comfort and acceptance of these masks.
- › Clear, soft vinyl for patient comfort and visual assessment.

CHIRO-OXY HIGH FLOW

Oxygen Mask High Flow (Adult/Pediatric/Neonatal)



- › Low resistance check valve. Prevents patients from rebreathing gas already exhaled
- › Odorless vinyl body. Soft mask is latex free, comfortable for patients, and material is safely disposable.
- › Adjustable nose clip. Allows mask to be formed to patient for tighter seal.

CHIRO-NEB


Nebulizer Mask (Adult/Pediatric/Neonatal)



- › Clear, soft vinyl for patient comfort and visual patient assessment.
- › Latex-free
- › Adjustable nose clip assures comfortable fit.
- › Specifically design for aerosol therapy.

CHIRO-PRONG

Nasal Prongs (Adult/Pediatric/Neonatal)



- › Delivers safe, low/moderate concentrations of oxygen for extended periods of time.
- › The soft nasal part of the anatomical shape reduces the risk of damage to the mucosa.
- › Comes with kink-resistant supply tubing.

HME FILTER



- › HME filter is used on patient breathing circuits to avoid any cross-contamination.
- › Appropriate for spontaneous respiration patients with a tracheostomy tube.
- › Low resistance to flow reduces load on breathing.

CHIRO-CISER

Lung Exerciser 3 Balls



- › Wide inspiratory flow rate range from 600 – 1200cc/sec to suit the widest range of patients.
- › Removable inspiratory limb for cleaning convenience.
- › Compact design & break-resistant plastic.

CATHETER MOUNT



- › Catheter Mount is a medical accessory that is attached with the tracheal tube to the end of the anesthetic breathing system.
- › Catheter Mount attached with breathing circuit minimizes the jerk of breathing to ET or TT reducing trauma to the trachea.

BAIN CIRCUIT



- › The process includes patient inspiring fresh gas from the outer reservoir tube and expiring into reservoir tube.
- › Convenient, lightweight anesthesia delivery system.

› **Unit Contents :**

- 22mm Cap
- 22mm corrugated tube
- 22mm expiratory valve
- 2 litres antistatic rebreathing bag
- Standard connectors
- Green Tube



VENTILATOR CIRCUIT




**Ventilator Circuit
without Water Trap**



**Chiro-Venti Kit
(Ventilator Circuit + HME Filter
+ Catheter Mount)**



VENTILATOR CIRCUIT



**Ventilator Circuit
with Water Trap**



**Ventilator Circuit
with two Water Traps**





SURGERY & DRAINAGE

Surgical drains are catheters placed near the surgical wound to remove blood, pus, or other fluid and prevent it from accumulating in the body.

CHIRO-ADK



- > Specially designed kit for post-operative abdominal drainage.
- > Drainage catheter along with a collection bag.
- > Catheter with atraumatic eyes for efficient drainage.



UMBILICAL CORD CLAMP



- > Suitable for clamping the umbilical cord of new born baby immediately after the birth.



CHEST DRAIN CATHETER



With Trocar




Without Trocar



- › Suitable for effective drainage after cardio-thoracic & thoracic surgery.
- › With radio-opaque line and markings at every 2 cm from the last eye.

YANKAUER SUCTION SET




- › Moulded Yankauer suction handle with flexible kink resistant tube.
- › Soft flexible adaptors at both ends of the tube provide safe grip to the handle as well as suction source.




CHIRO-WOUND

Wound Drainage Set



Curved Needle with Matching Catheter



Spare Perforated Catheter



Bellow Unit Capacity 800 ml



Connecting Tube with Clamp and "Y" connector




chironsurgimed.com

- Close wound drainage device under negative pressure.
- Provided with the option to use one or two catheters post-operatively.
- Multiperforated catheters provided with radio opaque line and eyes for trauma free performance.

Sizes Available :



CAUTERY PENCIL




- > Cautery pencils are widely employed in surgical procedures for their ability to achieve precise and controlled cauterization.
- > 3 Metre long wire.



CHIRO-SHEET

Corrugated Drainage Sheet



- > Designed for efficient multichannel wound drainage specially where airtight closure of wound is not possible.






GASTROENTEROLOGY

Gastroenterology is the branch of medicine focused on the digestive system and its disorders.

RYLES TUBE




- › Designed for introduction of nutrition and aspiration of gastric contents.
- › Four lateral eyes are provided for efficient aspiration and administration.

Sizes Available :



FEEDING TUBE



- › Infant feeding tube with graduation.
- › Made from DEHP Free, non-toxic PVC material completely non-irritant to delicate mucosa.
- › The tube features low friction, super smooth frozen surface ensuring non-traumatic intubation.

Sizes Available :




MUCUS EXTRACTOR




- › Suitable for obtaining mucus specimen for microbiological examination.
- › Clear and transparent container allowing visual examination of mucus specimen.
- › Universal plain funnel connector for easy connection.

Sizes Available :



MISCELLANEOUS



Cannula Fixture



Microporous Tape



NOTES


chironsurgimed.com





CHIRON

CHIRON SURGIMED LLP

 883/1, GIDC Road No.3, Mansa - 382845,
Gujarat, India.

 www.chironsurgimed.com

 info@chironsurgimed.com

 +91 8758888150, +91 2763 272480

Designed By



+91 9560827828
adhar@studioad.works