



# PRODUCT CATALOGUE



# ABOUT US



**Chiron Surgimed Private Limited** is a pioneering company at the forefront of innovative medical technology. Established in 2013, we have been dedicated to revolutionizing healthcare by developing cutting-edge medical devices that enhance patient care and improve clinical outcomes.

Our company specializes in the design, manufacturing, and distribution of a wide range of medical devices catering to various medical specialties such as Infusion Therapy, Respiratory Care, Urology and more.

In essence, **Chiron Surgimed Private Limited** is not just a company; it's a commitment to improving lives through innovative medical technology.

# MEET THE TOP BRASS



## Mr. Naresh Patel

**Founder & Chairman**

Founder to several companies with background in speciality polymer and in-depth experience of over 3 decades. His core competency is plastic polymers in the field of medical devices. Extremely diligent and hardworking, he is highly devoted spiritually which reflects in his day-to-day activities. Believes in working hard and simple living.



## Mr. Abhay Kumar

**CEO & Managing Director**

He possesses strong leadership qualities with over 35 years of professional experience in leading global medical device companies. He is a post graduate in Microbiology and done business management from premiere institute. He has strong knowledge of global business development, product development, business management and quality management skills.



## Mr. Harsh Patel

**Director**

A young entrepreneur with zeal to excel in Business Development and Marketing. Graduate in Commerce and keenly involved in Business Development activities and Production. Extremely hard working and intelligent person who wants to take Chiron to the next level. Talented and passionate about learning all the tricks of the trade.

# CONTENTS

## INFUSION THERAPY

<b>CHIRO-DRIP OPEL</b>	<b>06</b>
<b>CHIRO-DRIP ASTRA</b>	<b>07</b>
<b>CHIRO-DRIP ACCUFUSION</b>	<b>08</b>
<b>CHIRO-DRIP VENTED</b>	<b>09</b>
<b>CHIRO-DRIP NON-VENTED</b>	<b>10</b>
<b>CHIRO-DRIP Y</b>	<b>11</b>
<b>CHIRO-DRIP ONCO</b>	<b>12</b>
<b>CHIRO-DRIP ONCO GUARD</b>	<b>13</b>
<b>CHIRO-GUARD</b>	<b>14</b>
<b>CHIRO-DRIP MICRO</b>	<b>15</b>
<b>CHIRO-DRIP MICRO Y</b>	<b>16</b>
<b>CHIRO-FLOW</b>	<b>17</b>
<b>CHIRO-BT</b>	<b>18</b>
<b>CHIRO-VOL</b>	<b>19</b>
<b>CHIRO-VOL Y</b>	<b>20</b>
<b>CHIRO-VEIN AND NEO</b>	<b>21</b>
<b>CHIRO-VEIN SAFETY</b>	<b>22</b>
<b>CHIRO-VEIN NEO SMALL WINGS</b>	<b>23</b>
<b>CHIRO-WAY</b>	<b>24</b>
<b>CHIRO-WAY LINE</b>	<b>24</b>
<b>CHIRO-ACCUSAFE LINE</b>	<b>25</b>
<b>CHIRO-XLINE</b>	<b>26</b>
<b>CHIRO-PML</b>	<b>27</b>

## UROLOGY & NEPHROLOGY

<b>CHIRO-CATH</b>	<b>29</b>
<b>CHIRO-BAG</b>	<b>30</b>
<b>CHIRO-METER</b>	<b>31</b>
<b>NELATON</b>	<b>32</b>
<b>URON</b>	<b>32</b>
<b>DYRON</b>	<b>33</b>
<b>CHIRO-CVC KIT</b>	<b>34</b>

## RESPIRATORY CARE

<b>CHIRO-AMBU</b>	<b>36</b>
-------------------	-----------

<b>CHIRO-GUEDEL</b>	<b>36</b>
<b>CHIRO-OXY</b>	<b>37</b>
<b>CHIRO-OXY HIGH FLOW</b>	<b>37</b>
<b>CHIRO-NEB</b>	<b>38</b>
<b>CHIRO-PRONG</b>	<b>38</b>
<b>CHIRO-HME</b>	<b>39</b>
<b>CHIRO-CISER</b>	<b>39</b>
<b>CHIRO-ETT UNCUFFED</b>	<b>40</b>
<b>CHIRO-ETT CUFFED</b>	<b>40</b>
<b>CHIRO-SUCTION</b>	<b>41</b>
<b>CHIRO-SUCTION THUMB CONTROL</b>	<b>41</b>
<b>CHIRO-MOUNT</b>	<b>42</b>
<b>CHIRO-BAIN</b>	<b>42</b>
<b>CHIRO-VENTI WITHOUT TRAP</b>	<b>43</b>
<b>CHIRO-VENTI KIT</b>	<b>43</b>
<b>CHIRO-VENTI WITH WATER TRAP</b>	<b>44</b>
<b>CHIRO-VENTI WITH DOUBLE WATER TRAPS</b>	<b>44</b>

## **SURGERY & DRAINAGE**

<b>CHIRO-ADK</b>	<b>46</b>
<b>CHIRO-CLAMP</b>	<b>46</b>
<b>CHIRO-CATH DRAINAGE</b>	<b>47</b>
<b>CHIRO-YS</b>	<b>47</b>
<b>CHIRO-WOUND</b>	<b>48</b>
<b>CHIRO-PENCIL</b>	<b>49</b>
<b>CHIRO-SHEET</b>	<b>49</b>

## **GASTROENTEROLOGY**

<b>CHIRO-RYLES</b>	<b>51</b>
<b>CHIRO-FEED</b>	<b>51</b>
<b>CHIRO-EXTRACTOR</b>	<b>52</b>

## **MISCELLANEOUS**

<b>CHIRO-FIX</b>	<b>52</b>
<b>CHIRO-PORE</b>	<b>52</b>
<b>CHIRO-CLENZ</b>	<b>52</b>



# INFUSION THERAPY

Infusion therapy is when you receive medication through a needle or catheter, usually intravenously (IV).

# CHIRO-DRIP OPEL

Advanced IV Set Vented

**New Launch**




Largest Chamber in the Segment



Roller Clamp with Greater Control



High Weightage Tubing



- › Featuring the highest volume chamber in its segment, it ensures optimal fluid delivery.
- › The precision-engineered roller clamp offers unmatched control.
- › The high-weightage tubing enhances durability and reliability.



# CHIRO-DRIP ASTRA

Advanced IV Set Non-Vented

**New Launch**




**Largest Chamber in the Segment**



**Roller Clamp with Greater Control**



**High Weightage Tubing**



- › Featuring the highest volume chamber in its segment, it ensures optimal fluid delivery.
- › The precision-engineered roller clamp offers unmatched control.
- › The high-weightage tubing enhances durability and reliability.



# CHIRO-DRIP ACCUFUSION

IV Set with Safe Infusion and Auto-Stop Feature

**New Launch**

CIV 0210149



Auto Stop Filter



Roller Clamp



Priming Filter



- › The Auto Stop filter is expertly designed to block air from entering the IV line when the bottle runs empty.
- › This filter also ensures a steady fluid level, significantly lowering the risk of air entering the IV line.

# CHIRO-DRIP

CIV 0341299

## IV Set Vented

Air Vented Spike



Wingless Variant Also Available



Roller Clamp





chironsurgimed.com

- › Sharp air vented spike.
- › Suitable for pressure & gravity infusions.
- › Latex-free.
- › DEHP-free.

# CHIRO-DRIP

## IV Set Non-Vented



CIV 1341299



**Non-Winged**



**Winged**




CIV 1341199

# CHIRO-DRIP Y

CIV 0241012

## IV Set Vented with Y Injection Site

Wingless Variant Also Available



Air Vented Spike



Y Injection Site





- › Sharp air vented spike.
- › Suitable for pressure & gravity infusions.
- › Latex-free.
- › DEHP-free.
- › Y-injection port designed for multiple injections.


# CHIRO-ONCO

## Oncology Infusion Set

CIV 0241012 ONCO



Y Injection Site



Air Vented Spike




› Oncology Infusion Set are designed for the administration of chemotherapy.


# CHIRO-ONCO GUARD

**Oncology IV Set with 0.2 Micron Filter with Y Injection Site**  
with option of Needle Free Connector

**CIV 0246810 ONCO**



0.2 Micron Filter



Roller Clamp



Needle Free Connector




- › Sharp air vented spike
- › Suitable for pressure & gravity infusions
- › Latex free
- › DEHP free
- › Y-injection port designed for multiple injections
- › 0.2 Micron Inline Filter

# CHIRO-GUARD

CIV 0246810

**IV Set Vented with 0.2 Micron Filter with Y Injection Site**  
with option of Needle Free Connector



Roller Clamp




0.2 Micron Filter



Needle Free Connector




- › Sharp air vented spike
- › Suitable for pressure & gravity infusions
- › Latex free
- › DEHP free
- › Y-injection port designed for multiple injections
- › 0.2 Micron Inline Filter



# CHIRO-DRIP MICRO

CIV 1347129

## Micro IV Set Non- Vented



Micro Drip



Roller Clamp



Latex Bulb  
with Hypodermic Needle




- › Micro I.V. Infusion set for paediatric use
- › Approximately 60 drops/ml
- › Sharp Spike for easy insertion



# CHIRO-DRIP MICRO Y

CIV 0248109

## Micro IV Set Vented with Y Injection Site



Y Injection Site



Micro Drip



chironsurgimed.com



- › Micro I.V. Infusion set for paediatric use
- › Approximately 60 drops/ml
- › Y-injection port designed for multiple injections
- › Sharp and Wingless Spike for easy insertion

# CHIRO-FLOW

## Flow Regulator

Needle Free Connector



CIV 0257810



CIV 2510999





## Flow Regulator

- › Ergonomic dial design.
- › Compatible with standard sets.
- › Y-injection port designed for multiple injections.
- › DEHP-Free soft and kink resistant PVC tubing.

**Flow Regulator integrated with IV Set**

# CHIRO-BT

## Blood Transfusion Set



CBT 0239999



**Non-Vented**

CBT 0249999




**Vented**

# CHIRO-VOL

CMV 2349999 : 110ML  
CMV 2359999 : 150ML

## Measured Volume Burette Set



Air Vented Spike



Silicone Floater



- › Burette chamber of 110ml/150ml.
- › Micro drip with drop size of 60 drops per ml.
- › DEHP-Free soft and kink resistant PVC tubing.

# CHIRO-VOL Y

CMV 0149999 : 110ML  
CMV 0159999 : 150ML

## Measured Volume Burette Set with Y Injection Site



Y Injection Site



Needle Free Connector





- › Burette chamber of 110ml/150ml.
- › Micro drip with drop size of 60 drops per ml.
- › Y-injection port designed for multiple injections.
- › DEHP-Free soft and kink resistant PVC tubing.



# CHIRO-VEIN AND NEO

## IV Cannula with Wings and Injection Port



### TECHNICAL SPECIFICATION


Size	Catheter OD(mm)	Catheter length(mm)	Flow Rate (Water) (ml/min)	Box Qty. (Inner/Outer)
14G	2.1	45	240	100/1000
16G	1.7	45	200	100/1000
17G	1.5	45	125	100/1000
18G	1.3	45	90	100/1000
20G	1.1	32	60	100/1000
22G	0.9	25	36	100/1000
24G	0.7	19	23	100/1000
26G	0.6	19	15	100/1000

**CCN 123XX99**

XX represents size from 16G to 26G

# CHIRO-VEIN SAFETY

## Safety IV Cannula with Wings and Injection Port




### TECHNICAL SPECIFICATION

Size	Catheter OD(mm)	Catheter length(mm)	Flow Rate (Water) (m/min)	Box Qty. (Inner/Outer)
18G	1.3	45	90	50/500
20G	1.1	32	60	50/500
22G	0.9	25	36	50/500
24G	0.7	19	23	50/500

**CCN 023XX99**  
XX represents size from 18G to 24G

# CHIRO-VEIN NEO SMALL WINGS

IV Cannula with Small Wings and without Injection Port



## TECHNICAL SPECIFICATION

Size	Catheter OD(mm)	Catheter length(mm)	Flow Rate (Water) (ml/min)	Box Qty. (Inner/Outer)
24G	0.7	19	15	100/1000
26G	0.6	19	15	100/1000


**CCN 12XX999**

XX represents sizes 24G and 26G

# CHIRO-WAY

CEL 0129999

## Three-Way Stopcock



CEL 012XXX9  
XXX REPRESENTS SIZE

## CHIRO-WAY LINE


### Three-Way Stopcock with Extension Tube

Length available from  
10cm to 300cm





### TECHNICAL SPECIFICATION

Options	Packaging Specification		
Tap Colour	Material	Tap Rotation	inner
Red	Lipid   Lipid Resistant	Smooth   45° Click	50
White	Lipid   Lipid Resistant	Smooth   45° Click	outer
Blue	Lipid   Lipid Resistant	Smooth   45° Click	1000



# CHIRO-ACCUSAFE LINE




chironsurgimed.com

## Multi-Infusion Extension Tube with Needle Free



# CHIRO-XLINE




**Extension Tube  
with 0.2 Micron Filter**

**Standard Extension  
Tube**

CEL 0120209

CEL 0120109



**Spiral Extension Tube** CEL 0123009

Length available from  
10cm to 300cm




# CHIRO-PML

CPM 012XXX9

XXX REPRESENTS SIZE

## Pressure Monitoring Line



Length Available from  
100cm to 200cm  
For Pressure Upto : 580 PSI/40 Bar



# UROLOGY & NEPHROLOGY


A nephrologist deals with problems of the kidney with medicines and a urologist deals with kidney ailments that require surgery. Any ailments in the urinary system, prostate, bladder and penis which require surgical intervention are tended to by a urologist.

# CHIRO-CATH

## Foley Catheter



### Latex Foley Catheter



Sizes Available :

■	8FR	CFL 0240899
■	10FR	CFL 0241099
■	12FR	CFL 0241299
■	14FR	CFL 0241499
■	16FR	CFL 0241699
■	18FR	CFL 0241899
■	20FR	CFL 0242099
■	22FR	CFL 0242299

Three-Way Variant also available.



### Silicon Foley Catheter




Sizes Available :

■	6FR	CFS 1240699
■	8FR	CFS 1240899
■	10FR	CFS 1241099
■	12FR	CFS 1241299
■	14FR	CFS 1241499
■	16FR	CFS 1241699
■	18FR	CFS 1241899

# CHIRO-BAG

## Urine Bag



CUB 0299999

Premium version features a high-quality plastic handle

CUB 1349999

Bottom drain variant introducing soon !  
CUB 1346999


- › Urine collecting bag used for urine drainage.
- › Two litre capacity bag provided with super smooth kink resistant tubing.
- › Efficient non-return valve with top outlet.



# CHIRO-METER

CUB 1345999

## Urometer



chironsurgimed.com



- › Graduated measured volume meter for accurate measurement of the urine output.
- › Urometer is directly attached to a 2 litre capacity urine bag having overflow facility.
- › Soft and kink resistant PVC tubing.



# NELATON

## Urethral Tube



**Sizes Available :**  
**from 8FR to 16FR**

**CNL OXX9999**  
XX REPRESENTS SIZE  
FROM 8 TO 16



- › Designed for short term bladder catheterization through urethra.
- › Manufactured from medical grade PVC compatible with catheter lubricants.




# URON

## DJ Stent



**OPEN END : CDJ OXX9999**  
**CLOSE END : CDJ 1XX9999**  
XX REPRESENTS SIZE  
FROM 5FR TO 7FR



- › Short and Long Term connection of Urethra and Bladder.
- › Closed Tip and Open Tip

# DYRON

CHD 1299999: DOUBLE LUMEN CATHETER  
CHD 1399999: TRIPLE LUMEN CATHETER

## Long-term and Short-term Haemodialysis Double/Triple Lumen Catheter




› Haemodialysis catheter are double/ triple lumen catheter that provide temporary vascular access for haemodialysis untill a permanent access is available or until another type of dialysis therapy is substituted.

### › Kit Contents :

Indwelling catheter  
Y-Introducer needle with valve  
J-tip guide wire  
Dilators- 2 Nos  
Luer lock syringe  
Scalpel

# CHIRO-CVC KIT

## Central Veinous Catheter Kit in Double/Triple Lumen Variant



› For administration of hyperosmolar solutions, measuring CVP, haemodialysis, plasmapheresis, rapid infusion of fluid among other uses.

### › Kit Contents :

- Indwelling catheter
- Introducer needle
- Guide wire with Advancer
- Syringe
- Syringe Guide
- Injection Caps Heparin
- Scalpel
- Dilator
- Patient Job Card



# RESPIRATORY CARE

Respiratory therapy is a field of treatment that helps individuals optimize their respiratory function, breathe more easily, and live more comfortably.

# CHIRO-AMBU

## Silicone Resuscitator (Adult/Pediatric/Neonatal)



CAB 0199999 : Adult  
CAB 0299999 : Pediatric  
CAB 0399999 : Neonatal



- > Silicone resuscitator (except mask, oxygen tubing and reservoir bag) can be autoclaved repeatedly.
- > It is semi-transparent and come with a pressure limitation valve for patient safety.

# CHIRO-GUEDEL

CGA XXX9999  
XXX REPRESENTS SIZE AVAILABLE

## Guedel Airways



- > The center passage permits the passage of suction catheters to clear secretions from the back of the throat.




Sizes Available :



000 00 0 1 2 3 4 5

# CHIRO-OXY

## Oxygen Mask (Adult/Pediatric)




COX 0299999 : Adult  
COX 1399999 : Pediatric

- › The clinically clean oxygen masks are used to administer oxygen.
- › The adjustable nose clip and the under-chin design increase the comfort and acceptance of these masks.
- › Clear, soft vinyl for patient comfort and visual assessment.

# CHIRO-OXY HIGH FLOW

## Oxygen Mask High Flow (Adult/Pediatric)



COX 0245999 : Adult  
COX 1345999: Pediatric

- › Low resistance check valve. Prevents patients from rebreathing gas already exhaled
- › Odorless vinyl body. Soft mask is latex free, comfortable for patients, and material is safely disposable.
- › Adjustable nose clip. Allows mask to be formed according to patient for tighter seal.

# CHIRO-NEB

CNM 0249999 : Adult  
CNM 1349999: Pediatric

## Nebulizer Mask (Adult/Pediatric)



**T-Piece Connector**




- › Clear, soft vinyl for patient comfort and visual patient assessment.
- › Latex-free
- › Adjustable nose clip assures comfortable fit.
- › Specifically design for aerosol therapy.



# CHIRO-PRONG

## Nasal Prongs (Adult/Pediatric/Neonatal)



CNC 0349999 : Adult  
CNC 1349999 : Pediatric  
CNC 2499999 : Neonatal



- › Delivers safe, low/moderate concentrations of oxygen for extended periods of time.
- › The soft nasal part of the anatomical shape reduces the risk of damage to the mucosa.
- › Comes with kink-resistant supply tubing.



# CHIRO-HME

CHF 1239999

## HME Filter



- › HME filter is used on patient breathing circuits to avoid any cross-contamination.
- › Appropriate for spontaneous respiration patients with a tracheostomy tube.
- › Low resistance to flow reduces load on breathing.

# CHIRO-CISER

CCC 0129999

## Lung Exerciser 3 Balls




- › Wide inspiratory flow rate range from 600 – 1200cc/sec to suit the widest range of patients.
- › Removable inspiratory tube for cleaning convenience.
- › Compact design & break-resistant plastic.

# CHIRO-ETT

**CET OXX9999**  
XX REPRESENTS SIZE

## Endotracheal Tubes Uncuffed

New Launch



- › Endotracheal tubes are medical devices used to maintain an open airway and facilitate breathing in patients who are unable to do so on their own.

**Sizes Available :**


From 2mm ID to 5mm ID

# CHIRO-ETT

**CET 1XX9999**  
XX REPRESENTS SIZE

## Endotracheal Tubes Cuffed

New Launch



- › Cuffed endotracheal tubes are commonly used in adults and sometimes in pediatric patients, especially in situations where long-term ventilation is needed or where there is a high risk of aspiration.

**Sizes Available :**


From 7mm ID to 10mm ID

# CHIRO-SUCTION

**CSC 0XX9999**

XX REPRESENTS SIZE FROM 8FR TO 18FR

## Suction Catheter



- › Suction Catheters are suitable for removal of secretion from mouth, trachea and bronchial tubes.
- › Provided with universal funnel shape connector for safe connection to standard suction equipment.

### Sizes Available :




# CHIRO-SUCTION

**CSC 1XX9999**

XX REPRESENTS SIZE FROM 8FR TO 18FR

## Suction Catheter with Thumb Control Connector



- › Suitable for power controlled extraction of bodily secretions, such as mucus from the upper airway.

### Sizes Available :



# CHIRO-MOUNT

CCM 0999999

## Catheter Mount



- › Catheter Mount is a medical accessory that is attached with the tracheal tube to the end of the anesthetic breathing system.
- › Catheter Mount attached with breathing circuit minimizes the jerk of breathing to ET or TT reducing trauma to the trachea.

# CHIRO-BAIN

## Bain Circuit



CBC 0999999 : Adult  
CBC 1999999 : Pediatric

- › The process includes patient inspiring fresh gas from the outer reservoir tube and expiring into reservoir tube.
- › Convenient, lightweight anesthesia delivery system.

### › Unit Contents :

22mm Cap  
22mm corrugated tube  
22mm expiratory valve  
2 litres antistatic rebreathing bag  
Standard connectors  
Green Tube

# CHIRO-VENTI


## Ventilator Circuit

CVC 0149999 : Adult  
CVC 0159999 : Pediatric



## Ventilator Circuit without Water Trap


CVK 0129999 : Adult  
CVK 0139999 : Pediatric




## Chiro-Venti Kit (Ventilator Circuit + HME Filter + Catheter Mount)

# CHIRO-VENTI

## Ventilator Circuit



CVC 0249999: Adult  
CVC 0259999 : Pediatric



## Ventilator Circuit with Water Trap

CVC 0349999: Adult  
CVC 0359999 : Pediatric  
CVC 0369999 : Neonatal



## Ventilator Circuit with two Water Traps





# **SURGERY & DRAINAGE**

**Surgical drains are catheters placed near the surgical wound to remove blood, pus, or other fluid and prevent it from accumulating in the body.**

# CHIRO-ADK

**CAK 012XX99**  
XX REPRESENTS SIZE  
FROM 20FR TO 32FR

## Abdominal Drainage Kit



- › Specially designed kit for post-operative abdominal drainage.
- › Drainage catheter along with a collection bag.
- › Catheter with atraumatic eyes for efficient drainage.



# CHIRO-CLAMP

**CUC 0999999**

## Umbilical Cord Clamp



- › Suitable for clamping the umbilical cord of new born baby immediately after the birth.



# CHIRO-CATH DRAINAGE

## Chest Catheter Drainage



**With Trocar**

**CCD OXX9999**  
XX REPRESENTS SIZE  
FROM 20FR TO 32FR



**Without Trocar**

**CCD 1XX9999**  
XX REPRESENTS SIZE  
FROM 20FR TO 32FR



- › Suitable for effective drainage after cardio-thoracic & thoracic surgery.
- › With radio-opaque line and markings at every 2 cm from the last eye.

# CHIRO-YS

**CYS-0239999**


## Yankauer Suction Set




- › Moulded Yankauer suction handle with flexible kink resistant tube.
- › Soft flexible adaptors at both ends of the tube provide safe grip to the handle as well as suction source.

# CHIRO-WOUND

## Wound Drainage Set



Curved Needle with  
Matching Catheter



Spare Perforated Catheter



Bellow Unit  
Capacity 800 ml



Connecting Tube with Clamp  
and "Y" connector



- › Close wound drainage device under negative pressure.
- › Provided with the option to use one or two catheters post-operatively.
- › Multiperforated catheters provided with radio opaque line and eyes for trauma free performance.



### Sizes Available :

				
10FG	12FG	14FG	16FG	18FG
CCW 0161099	CCW 0271299	CCW 0381499	CCW 0410169	CCW 0511189

# CHIRO-PENCIL

CCP 0999999

## Cautery Pencil






- › Cautery pencils are widely employed in surgical procedures for their ability to achieve precise and controlled cauterization.
- › 3 Metre long wire.



# CHIRO-SHEET

CDS 0999999

## Corrugated Drainage Sheet



- › Designed for efficient multichannel wound drainage specially where airtight closure of wound is not possible.





# GASTROENTEROLOGY


**Gastroenterology is the branch of medicine focused on the digestive system and its disorders.**

# CHIRO-RYLES

**CRT OXX9999**

XX REPRESENTS SIZE FROM 8FR TO 18FR

## Ryle's Tube



- › Designed for introduction of nutrition and aspiration of gastric contents.
- › Four lateral eyes are provided for efficient aspiration and administration.

### Sizes Available :




# CHIRO-FEED

**CFT OXX9999**

XX REPRESENTS SIZE FROM 5 FG TO 10 FG

## Infant Feeding Tube



- › Infant feeding tube with graduation.
- › Made from DEHP Free, non-toxic PVC material completely non-irritant to delicate mucosa.
- › The tube features low friction, super smooth frozen surface ensuring non-traumatic intubation.

### Sizes Available :



# CHIRO-EXTRACTOR

CME 012XX99  
XX REPRESENTS SIZE FROM 6FR TO 18FR

## Mucus Extractor



- › Suitable for obtaining mucus specimen for microbiological examination.
- › Clear and transparent container allowing visual examination of mucus specimen.
- › Universal plain funnel connector for easy connection.

### Sizes Available :



6FR 8FR 10FR 12FR 14FR 16FR 18FR



## MISCELLANEOUS



CCF 09999999



**CHIRO-FIX**  
Cannula Fixture



**CHIRO-PORE**  
Microporous Tape

CST 00599999 : 5 METER  
CST 00999999 : 9 METER



**CHIRO-CLENZ**  
Bed Bath Wipes



# NOTES

[illegible]

# NOTES

[illegible]

# NOTES

[illegible]

# NOTES

[illegible]





# CHIRON

## CHIRON SURGIMED PRIVATE LIMITED

 Survey No. - 883/1, G.I.D.C Road No. - 3,  
Mansa - 382845, Dist- Gandhinagar, Gujarat (India)

 [www.chironsurgimed.com](http://www.chironsurgimed.com)

 [enquiry@chironsurgimed.com](mailto:enquiry@chironsurgimed.com)

 +91 6357000946



Scan to Download  
Catalogue

Designed By  
**ADHAR  
GUPTA  
STUDIO**  
+91 9560827828  
[me@adhargupta.com](mailto:me@adhargupta.com)